

DESIGNED FOR ULTRASOUND RESEARCH



- The World's First Ultra High to Low Frequency Imaging System \lhd
 - Open and Configurable Architecture <
 - High Resolution Real-Time Imaging \lhd

FUJIFILM | VISUALSONICS

Vevo F2

Introducing the Vevo F2

The World's First Ultra High to Low Frequency Imaging System for Ultrasound Research

For over twenty years, FUJIFILM VisualSonics has been delivering the best-in-class, ultra high frequency ultrasound and photoacoustic imaging solutions to the scientific research community. With the Vevo F2, we now expand our reach to satisfy the imaging needs of acoustic researchers, ultrasound engineers and those that may benefit from both ultra high to low frequency imaging capabilities.



Flexible

Ultra high to low frequency imaging (71-1 MHz)

1	~	
	3	
	2	
	TATA	

One System Adaptable for imaging small to large animals



Intuitive

Easy-to-use graphical interface

individual channel data (VADA)

Open Architecture

Access pre-beamformed

High Resolution Real-Time Imaging



Power Doppler Mode

Visualize and quantify blood flow



3D-Mode

Accurately quantify volumes-can be combined with Power Doppler or photoacoustics

Imagine the possibilities:

- Plane Wave Imaging
- Ultrafast Doppler
- Signal Processing and Beamforming

Small to Large Animal Imaging

▶ Tissue Characterization

Super Resolution Imaging

Photoacoustics

Compatible with the Vevo LAZR-X laser cart for multi-modal imaging

The Vevo F2 Imaging Platform



ADC resolution up to 12 bits

4

Open and Configurable Architecture

Vevo Advanced Data Acquisition (VADA)

Equipped with Vevo Advanced Data Acquisition (VADA), the Vevo F2 allows access to pre-beamformed individual channel data via an all new, easy-to-use, graphical interface. With full control over transmit profiles, researchers now have the power and freedom to develop and explore new imaging methods in a quick, iterative fashion-going beyond existing imaging modes.



Preview

-

3



Pulse Sequence



Transmit Delay



Configurable by User VADA (Vevo Advanced Data Acquisition)



Multi-Channel Data

World's First High to Low **Frequency Imaging System**

The Vevo F2 offers an expanded range of frequencies (71-1 MHz). Users now have the flexibility to image at low frequency for penetration and ultra high frequency for resolution using one platform.







UHF22x 22-10 MHz

L38xp 10-5 MHz





UHF57x 57-25 MHz



molecular imaging, co-registered with high

Target Areas of Research



Plane-wave Implementation

Implement plane-wave techniques for ultrafast ultrasound imaging for applications such as ultrafast Doppler and super-resolution ultrasound



Beamforming Algorithm Development

Test novel beamforming techniques for image reconstruction



External Devices Syncing

Coordinate timing between HIFU pulses for imaging, or shear wave generation for elastography measurements



UHF71x

71-30 MHz

Small to Large Animals

Conduct imaging and analysis of small and large animals on one platform to streamline data collection

Vevo F2

Vevo F2 Transducers

High to Low: Flexibility at your Fingertips

The Vevo F2 Imaging System is compatible with a greater range of transducers than ever before. With this new expanded range (71-1 MHz), users now have the flexibility to image a broad spectrum of animals from small to large, using the same imaging platform!

Transition from one transducer to another quickly; the Vevo F2 allows for **three transducers** to be connected with one transducer active at any given time.

High Frequency Transducers

Model	Туре	Bandwidth	Possible uses <i>in vivo</i>
			Mouse embryology
UHF71x	Linear	71-30 MHz	Vascular and epidermal imaging
			Ophthalmology
UHF57x	Linear	57-25 MHz	Mouse cardiovascular, abdominal, reproductive
			Mouse/rat embryology
			Small rat vascular
UHF46x	Linear	46-20 MHz	Mouse cardiovascular
			Rat abdominal
			Rabbit ophthalmology
			Rat/rabbit vascular
UHF29x	Linear	29-15 MHz	▷ Rat cardiology and abdominal (<250 g)
UHF22x	Linear	22-10 MHz	\triangleright Rat cardiology and abdominal (<500 g)
			Rabbit cardiovascular

Low Frequency Transducers

Model	Туре	Bandwidth	Possible uses in vivo		
L38xp	Linear	10-5 MHz	Large animal abdominal		
			Rat and rabbit cardiovascular		Proudly pa
			Low frequency photoacoustic imaging	Ь	FUJIFILM S
P10xp	Phased	8-4 MHz	Large animal cardiology		offer low fi
			Large animal abdominal		transducer

rtnering with onosite to equency

-1

High Quality Imaging for Precise Visualization of Tissue Structures



Mouse carotid artery UHF71x (71-30 MHz)



Mouse Spleen & Kidney Vasculature UHF57x (57-25 MHz)





Mouse Brain Vasculature UHF29x (29-15 MHz)



Diseased rat liver UHF22x (22-10 MHz)

Vevo F2

10

Vevo LAB Analysis Software

Data management and analysis with Vevo LAB workstation software. Export data to other third-party data processing tools.



Vevo Imaging Station

Standardize image acquisition and quantification to ensure repeatable, reproducible results and high-throughput workflow for multiple animal studies.

- ▶ Warmed platform for maintaining optimal physiological conditions for small animals
- Integrated & displayed physiological monitoring: ECG, heart rate, core temperature, respiration and blood pressure
- ▷ Transducer mounting system for precision and hands-free scanning
- Precision micro-injection system for injections or extraction procedures
- Compatible with the Vevo Compact Anesthesia System, the Vevo Whole Body Imaging Setup, the Vevo E-Box and Vevo BRAIN

Accessories



Anesthesia System



Vevo BRAIN Stereotactic Frame & Atlas



Vevo Whole Body Imaging Setup

Vevo Technology Timeline



Vevo Support

The advanced technology of the Vevo F2 high resolution imaging platform is accompanied by an integrated approach to service and support.

Applications support and training customized to your needs

On-site customer training

- Customized hands-on
- education

Online resources

▶ Live & on-demand webinars

- Imaging guides
- Video tutorials
- Grant support program
- Publications libraries
- Image galleries
- Exclusive customer
- resource portal

For additional resources, support or service requests, visit our website: visualsonics.com

- Technical and scientific support

On-site and online

- support
- Scientific application expertise

"The possibility of working with low and high frequency is really interesting. The VADA interface is intuitive and very easy to work with. The combination of VADA and the broadband frequencies are clear differentiators of Vevo F2 over other providers."

Magnus Cinthio, Associate Professor

Lund University

FUJIFILM VISUALSONICS

Seeing More Matters visualsonics.com

VisualSonics, the VisualSonics logo, Vevo, Vevo MicroMarker and EKV are registered and unregistered trademarks of FUJIFILM Sonosite, Inc. in various jurisdictions.

FUJIFILM is a trademark and registered trademark of FUJIFILM Corporation in various jurisdictions.

© 2021 VisualSonics Inc. All rights reserved. FUJIFILM Visualsonics, Inc is a subsidiary of FUJIFILM Sonosite, Inc.

MKT03495 (Rev 1.1)