

ARCHER VARIANTPlex™ COMPLETE SOLID TUMOR AND PAN SOLID TUMOR

Easily identify key oncology biomarkers

The expansive landscape of genetic biomarkers relevant to oncology continues to grow due to diligent research efforts across the globe. Informed decisions not only rely on identifying oncology genes, but understanding the magnitude of genome-wide signatures such as:

- **Microsatellite instability (MSI)**, which measures changes in the number of repeats found in microsatellites that could indicate defects in DNA mismatch repair (MMR) genes
- **Tumor mutational burden (TMB)**, which represents the number of acquired mutations per megabase (mut/Mb) of sequenced DNA and functions as a predictive biomarker for neoantigen production by cancer cells.

Consequently, creating comprehensive genomic profiles (CGP) of cancers requires data across a broad genomic target space for all types of mutations. In addition, the data needs to flex with an ever-changing landscape of knowledge.

Archer VARIANTPlex panels for solid tumors provide comprehensive coverage, flexibility, and scalability using proven chemistry for fast answers to your questions.

VARIANTPlex panels for solid tumor research



Comprehensive

Identify relevant oncology genes and genomic signatures such as MSI and TMB



Flexible

Add custom content to your NGS panels for state-of-the-art assays*



Fast results

1.5-day library preparation with 3.5 hours hands-on-time



AMP™ chemistry

Proven chemistry for dependable NGS libraries



Archer™ Analysis

Powerful analysis platform for actionable insight

*Custom targets do not factor into TMB calculation due to algorithm constraints.

Solid tumor analysis for comprehensive or focused targets

VARIANTPlex Pan Solid Tumor


- CNVs
- SNVs
- Indels
- Microsatellite instability (MSI)
- Tumor Mutational Burden (TMB)[†]
- Covers coding region of 185 genes

VARIANTPlex Complete Solid Tumor


- CNVs
- SNVs
- Indels
- Microsatellite instability (MSI)
- Tumor Mutational Burden (TMB)[†]
- Covers coding region of 430 genes

[†]TMB analysis requires use of VARIANTPlex Pan Solid Tumor or VARIANTPlex Complete Solid Tumor. Custom TMB panels must contain all Pan or Complete Solid Tumor content.

A. Pan Solid Tumor



B. Complete Solid Tumor



	10 Mut/Mbp	(-) TMB-Intermediate 5–20 Mut/Mbp
Concordance	95.0%	100%
Positive % agreement (PPA)	93.2%	100%
Negative % agreement (NPA)	100%	100%

	10 Mut/Mbp	(-) TMB-Intermediate 6–17.5 Mut/Mbp
Concordance	83.8%	100%
Positive % agreement (PPA)	81.5%	100%
Negative % agreement (NPA)	100%	100%

Figure 1. VARIANTPlex TMB identification is highly concordant with orthogonal methods. 50–100 ng of DNA from samples with orthogonal TMB characterizations was converted into libraries using VARIANTPlex chemistry, sequenced on Illumina® NextSeq™ 500/550, NextSeq 2000, or NovaSeq® instruments, then processed through Archer Analysis. **(A)** Reference: 11 SeraCare® TMB standards and WES-characterized reference cell lines (ATCC or GIAB); tissue samples: 39 de-identified FFPE tumor tissue samples characterized with a 500+ gene hybrid capture NGS assay, 2 de-identified FFPE normal adjacent tissue samples, 4 de-identified peripheral blood samples. **(B)** Reference: 11 SeraCare TMB standards and 5 WES-characterized reference cell lines; tissue samples: 87 de-identified FFPE tumor tissue samples characterized with a 500+ gene hybrid capture NGS assay, 2 de-identified FFPE normal adjacent tissue samples, 10 de-identified peripheral blood samples. Normal cell lines and peripheral blood samples were assumed to have a TMB of 0. (-) TMB-Intermediate results exclude samples within the indicated default range for Archer Analysis across all assays evaluated.

Detect MSI with reliability

MSI is reliably identified with high concordance to orthogonal methods and is enabled for both VARIANTPlex Pan Solid Tumor and VARIANTPlex Complete Solid Tumor.

Panel coverage:

- 114 microsatellite loci targeted across the genome

Analytical Performance:

- 98.3% PPA (95% CI = 95.6–99.5%)
- 98.4% NPA (95% CI = 94.4–99.8%)

Analysis outputs:

- % unstable loci + MSS, intermediate, and MSI-H classifications




Figure 2. VARIANTPlex MSI identification is highly concordant with orthogonal methods. 50–250 ng of DNA from de-identified FFPE samples ($n = 362$) with orthogonal MSI characterizations was converted into libraries using VARIANTPlex chemistry and sequenced on Illumina NextSeq 500/550. Archer Analysis identified the percentage of unstable loci which is plotted on the y-axis; $>30\%$ = MSI-H (microsatellite instability-high); $20\text{--}30\%$ = intermediate; $<20\%$ = MSS (microsatellite stable). Orthogonal classification was determined by PCR fragment analysis or immunohistochemistry (IHC) to identify whether the sample was deficient in mismatch repair (dMMR) or MMR proficient.

Archer Analysis 7.2

Sensitive and accurate data analysis is essential for informed decisions. Archer Analysis has been updated to identify TMB and MSI signatures in addition to other variants in key genes associated with solid tumors.



More data in less time

MSI and TMB analysis

Faster analysis for all panels



Support for CGP

Comprehensive analysis

DNA and RNA in one platform



Architecture update

Improved analysis performance

Infrastructure updates provide more security and reliability

For a full list of gene targets, see archerdx.com/research-products/solid-tumor-research/

New expanded panels	Size	Recommended reads	Product #
VARIANTPlex Pan Solid Tumor	~669 kb	25 M	AB0142
VARIANTPlex Complete Solid Tumor	~1.42 Mb	45 M	AB0144

Don't see your desired targets?
You can design your own panel or add to an existing one at the Archer Assay Marketplace, just go to <https://assay-marketplace.archerdx.com/>
VARIANTPlex panels are customizable

Assay Marketplace

Easily identify key oncology biomarkers

The powerful combination of Archer VARIANTPlex panels and purpose-built Archer Analysis software unlocks new capabilities without compromising performance on established biomarkers.

For research use only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2023 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners, and may be registered in the USA and/or other jurisdictions. For specific trademark and licensing information, see www.idtdna.com/trademarks. Doc ID: RUO23-2109_001 07/23

