

## Your perfect partner for automatic DNA Flow Technology

Automatic DNA Flow hybridization, completely guided by hybriSoft software










- √ Fully automated
- √ UV light for DNA decontamination
- √ Extremely robust and compact
- √ Fast, sensitive and reliable diagnostic
- √ Thermostatic reaction chamber
- √ Vacuum pump control included
- √ Processed samples from 1 to 24 per run
- √ Execution time: 30-120min
- √ hybriSoft: specific software developed for hybriSpot platforms







REPORT
HH MOSITIVE
GENOTYPES DETECTED
High Risk
4
5
LOW Risk
7
7. 72

- ybriSpot platforms
  ▶ Samples administration
- Automatic capture of images
- Automatic analysis of final arrays
- Printed report

Reference: VIT-HS24

Developed and manufactured by



VITRO, S.A.
Luis Fuentes Bejarano, 60
Edificio Nudo Norte (Local 3)
41020 Sevilla
Telf: (+34) 95 493 32 00 / (+34) 91 382 16 20
export@vitro.bio
www.vitro.bio



