







md stainer

Automation for Immunohistochemistry
methods, FISH and CISH



Automatic

MD STAINER

Fully automated IHC, FISH and CISH

Capacity

36 determinations using independent and motorized incubation chambers. Allows overnight use, being able to process a total of 72 samples in one day

Quality

MD STAINER delivers high quality, reproducibility and consistency results

Identification

MD STAINER has an identification system by radiofrequency to control and verify the availability and status of each reagent. It is fully configurable by the user and provides information: batch, test number, etc.

Traceability

Better control and traceability of the sample.
A single labeling of the sample is performed, identifying patients through labels scan with 2D bar codes

Tracking

MD STAINER performs a variety of data tracking. Reports can be generated by patient information, reagent used or performed cycles

Security

MD STAINER works with low volume of waste, reduces the risk of exposure to hazardous chemicals and provides a safer working environment



Parque Tecnológico Ciencias de la Salud
Avda. del Conocimiento, 100
18016 Granada (Spain)
Tel: (+34) 958 27 14 49
Fax: (+34) 958 27 14 34
Hotline: (+34) 958 10 50 11
export@vitro.bio
www.vitro.bio



Distributed by