



AviPure™
Automated Extraction Solutions

Fully-Automated Nucleic Acid Extraction Solution



AviPure AX⁹⁶™
Automated Extraction Solutions



HIGH CAPACITY

- Dispensing volume: 20-1000 µl
- Real-time run monitoring



USER FRIENDLY OPERATION

- Easy to use Touch screen
- Built-in >42 programs
- Design - wide spacing silo configuration prevents cross-contamination



SMOOTH WORKFLOW

- Pre-aliquoted reagents - no direct contact with reagents/chemicals
- Sampling>Addition>Extraction



GREAT PERFORMANCE

- Rapid heating-cooling system minimizes waiting time and preserves sample integrity
- Stable, silent operation without vibration
- Automatic UV disinfection & HEPA filter system, avoid contamination





AviPure™
Automated Extraction Solutions

Fully-Automated Nucleic Acid Extraction Solution



AviPure Rx™
Automated Extraction Solutions



Automated RNA/DNA purification



Magnetic Bead-based extraction - efficient RNA/DNA binding



Pre-aliquoted reagents - eliminates contamination risk



Hands-on Time - maximum ~30-40min



Low sample requirement - just 200µl



Sample types incl. serum, plasma, naso-/oro-pharyngeal swabs

ORDER INFORMATION

REFERENCE NO.

DESCRIPTION

AVP-96

AviPure AX 96
Fully-Automated Nucleic Acid Extraction System

AVP10196

AviPure Rx
Nucleic Acid Extraction Reagents Kit

96 tests/Kit

